



## Invitron Total Proinsulin Assay Kit

IV2-003/103

### Introduction

The Invitron Total proinsulin Assay is an immunometric assay for the sum quantitation of intact proinsulin, 32-33 split proinsulin and des 31-32 split proinsulin in human plasma samples.

### Summary

Proinsulin is a precursor molecule for insulin and is synthesised by the pancreatic  $\beta$ -cells. Under normal circumstances, virtually all proinsulin is cleaved at residues 32-33 and 65-66 to produce insulin during the formation of secretory granules. Some proinsulin is released into the circulation, mainly as des 31-32 split proinsulin together but also with lesser amounts of 32-33 split proinsulin and intact proinsulin. It is generally believed that 65-66 split proinsulin and des 64-65 split proinsulin do not appear in detectable amounts in the circulation. Increased concentrations of circulating proinsulin may occur in non-insulin dependent (type II) diabetes and in patients with insulinoma. When used in conjunction with a specific intact proinsulin assay, measurements of total proinsulin can provide useful information on changes in the processing of insulin in such situations.

### Principle

The Invitron Total Proinsulin Assay is a two-site immunoassay, employing a specific solid phase antibody immobilised on microtitre wells, and a soluble antibody labelled with a chemiluminescent acridinium ester. The sample is incubated in the microtitre well together with a buffer and, after a wash step, the labelled antibody solution is added. This is followed by a further wash step to remove unbound labelled antibody before measurement. The bound luminescence is quantified by a microtitre plate luminometer capable of in situ reagent addition. The luminescent reaction is a rapid flash type (>95% complete in 1 second) which permits the entire plate to be read in approximately 5 minutes.

Rev 01

1

### Kit Contents

- **Standards:** Recombinant intact proinsulin in a serum matrix, lyophilised and sealed under vacuum for stability.
- **Labelled Antibody Concentrate:** Chemiluminescent labelled antibody in a protein matrix including preservatives and 0.05% sodium azide.
- **Labelled Antibody Diluent:** For diluting the labelled antibody to its working strength. Protein matrix including preservatives and 0.05% sodium azide.
- **Coated Microtitre Plate:** Microtitre plate coated with a specific monoclonal antibody. The plate is sealed inside a foil pouch with a dessicant to maintain a moisture-free environment.
- **Wash Buffer Concentrate:** (x10) phosphate buffered saline containing a detergent and 0.09% sodium azide.
- **Plate sealers**
- **Product Insert**


### Materials Required But Not Provided

- Deionised water
- Detection reagents (Invitron Catalogue No. IV1-001)
- Uncoated strips

### Procedure

1. Bring all kit components and samples to room temperature before use.
2. Reconstitute each of the standards by the addition of 1 ml of deionised water. Allow these to stand for 5 minutes, then mix gently to ensure all solids are dissolved.
3. Pipette 900  $\mu$ l of labelled antibody concentrate into one bottle of labelled antibody diluent and mix thoroughly.
4. Make up working strength wash buffer by diluting 1 part of wash buffer concentrate with 9 parts of deionised water.
5. Assemble the required number of coated strips in the plate holder. Any strips not used immediately may be stored inside a sealed polythene bag with silica gel dessicant. Make sure to fill remaining spaces in the plate holder with uncoated strips to ensure uniform heat transfer during incubation.
6. Pipette 50  $\mu$ l sample buffer into each well.
7. Pipette 50  $\mu$ l standard or sample into the wells. Standards must be run in duplicate.
8. Attach the plate sealer and incubate for 2 hours at 37°C.
9. Remove the plate sealer and perform 3 wash cycles with working strength wash buffer using an automatic plate washer.
10. Pipette 100  $\mu$ l labelled antibody solution into each well.
11. Attach the plate sealer and incubate for a further 1 h at 37°C.
12. Remove the plate sealer and perform 3 wash cycles with working strength wash buffer using an automatic plate washer.
13. Measure the light output from each well in a plate luminometer.

## Typical Standard Curve



N.B. This curve is for illustration purposes only, and must not be used for result calculation. RLU = Relative Light Units.

## Precision Profile

The precision of duplicate measurements was calculated for 86 patient samples. The mean coefficient of variation for these duplicates, which covered the analytical range 5-235 pmol/l was 6.1%.

## Between Assay Precision

| Total Proinsulin (pmol/l) | CV%  | n |
|---------------------------|------|---|
| 33.8                      | 12.8 | 7 |
| 81.9                      | 13.4 | 7 |
| 176                       | 6.2  | 7 |

## Spiking Recovery

Four plasma samples containing low endogenous proinsulin were spiked with 49.5 and 98.0 pmol/l recombinant proinsulin (samples 1 & 4) or with 98.0 and 192 pmol/l recombinant proinsulin (samples 2 & 3). Recoveries are shown as percentages of the expected result.

| Sample    | 1   | 2   | 3   | 4  |
|-----------|-----|-----|-----|----|
| 1st Spike | 102 | 104 | 102 | 94 |
| 2nd Spike | 99  | 98  | 91  | 92 |

Mean spiking recovery was 97.8%.

2

3

## Dilution Recovery

Five plasma samples containing 55-160 pmol/l proinsulin were diluted with zero standard by factors of 12.5%, 25% and 50%. Recoveries of proinsulin are shown as percentages of the expected result.

| Sample         | 1   | 2   | 3   | 4   | 5  |
|----------------|-----|-----|-----|-----|----|
| 12.5% dilution | 116 | 100 | 95  | 93  | 94 |
| 25% dilution   | 98  | 111 | 105 | 96  | 88 |
| 50% dilution   | 105 | 115 | 107 | 100 | 96 |

Mean dilution recovery was 101.3%.

## Cross Reactivity

Cross reactivity of related proteins were investigated at concentrations of 200 pM. Results are expressed as percentages of the reactivity of an identical concentration of intact proinsulin.

| Peptide                    | CR (%) |
|----------------------------|--------|
| Intact proinsulin          | 100    |
| 32-33 split proinsulin     | 97     |
| Des 31-32 split proinsulin | 100    |
| 65-66 split proinsulin     | 12     |
| Des 64-65 split proinsulin | 1.0    |
| Insulin                    | 2.2    |
| C-peptide                  | 0.0    |

## Sensitivity

Sensitivity was estimated as two standard deviations from the mean of 20 replicates of a zero standard. Calculated in this way, analytical sensitivity of the Total Proinsulin Assay is 0.35 pmol/l

## Range of Standards (Typically)

0-250 pmol/l

## Specimen

Plasma samples should be used for the total proinsulin assay.

## High Dose Hook Effect

Because of the assay architecture, which employs separate incubations with solid phase and labelled antibodies, no high dose hook effect is experienced.

4

5

## Safety Precautions

1. Use good laboratory practice at all times.
2. Wear appropriate protective clothing when using this kit and when handling samples.
3. This kit contains no human-derived material.

## Technical Precautions

1. Do not use components after their expiry date.
2. Once components have been opened or reconstituted, they can be used within a two-week period, provided they have been stored at 2-8°C.
3. Do not mix reagents from different lots.
4. Take care to avoid contamination of any reagent.

**This kit is intended for research use only.**

For additional information and product support please contact:

Invitron Limited  
Wyastone Business Park  
Wyastone Leys  
Monmouth NP25 3SR  
United Kingdom

Tel: +44 (0)1600 891536

Fax: +44 (0)1600 891537

Email: [info@invitron.co.uk](mailto:info@invitron.co.uk)

[www.invitron.co.uk](http://www.invitron.co.uk)



**Invitron Limited**  
Wyastone Business Park  
Wyastone Leys  
Monmouth NP25 3SR  
United Kingdom

Tel: +44 (0)1600 891536  
Fax: +44 (0)1600 891537  
Email: [info@invitron.co.uk](mailto:info@invitron.co.uk)  
[www.invitron.co.uk](http://www.invitron.co.uk)